



APEL CO., LTD. JAPAN

**APEL Hemoglobinmeter
MADE IN JAPAN**

HG-220

Hemoglobinmeter



Compact and simple

No lamp replacement

Thanks to LED, No Lamp replacement is required.

Auto-calibration

Thanks to Auto-calibration, it is easy to start HG-220 for daily measurement.

Compact & light weight

HG-220's foot print is small and weights around 0.6kg.





HG-220

Hemoglobinmeter

Specifications:

Measurement method	Cyanmethemoglobin method
Display	3-1/2 digit(7 Segment) LED
Wavelength	540 nm
Measuring Range	g/dl : 0.0 - 30.0 mmol/l : 0 - 18.6
Light source	Ultra-high brightness green LED
Photo detector	Silicon-photodiode
Sample volume	1.0 ml(minimum)
Accuracy	≤ ± 1.5 %
Dilution ratio	1 : 251
Power	AC Adapter (Input AC 100 - 240 VAC, Output DC12V 5A)
Frequency	50 - 60 Hz
Dimensions	117(W)x 185(D) x 77(H)mm
Cuvette	Square glass cuvette(12.5 x 12.5 x 45mm, outer size)
Weight	Net 0.6 kg (Main body)

Specifications and appearance are subject to change without notice.

Standard Accessories:

Operation manual	1 copy
AC adapter + cord (Apel's version)	1 pc
12.5 x 12.5 x 45 mm square cuvette	2 pcs

MADE IN JAPAN



APEL CO., LTD. JAPAN

Website

<http://www.apel.co.jp>

Distributor